



SPEED



POWER



POSTURE



SYMMETRY

APPLICATION ANDROID TABLET



BLUETOOTH CONNECTION



SCREEN MIRRORING





MULTICOLOR CHARTS



ZOOMABLE SCALE



INDEPENDENT SETUP OF MODULS

DYNAMIC VECTORING


FREEDOM OF MOVEMENT

Thanks to this patented technology, the load of exercise is always perpendicular to the ground or in other directions for any training requirements. It is possible to set a variable angle from 0° to 45°, mainteining then always constant the vector component, or set a fixed point of the cable exit. This allows you to fit the natural biomechanics of human body.










IN VERTICAL SPACE SAVER WORKBASE

The two tubolar modules can be **applied directly to a wall**, thanks to the **simple brackets** provided with rapid attachment/detachment.

You can pre-install additional brackets in other environments and use the same modules without limitations.


This configuration allows to easily optimize the space.



SOME HORIZONTAL WORKOUT EXAMPLES

















SOME VERTICAL WORKOUT EXAMPLES



IN HORIZONTAL A SOLID WORKBASE

Thanks to the platform accessory, it is possible to position the **tubolar modules horizontally** with the same **rapid** fastening/unfastening system.

In this configuration it is possible to perform multiple exsercises in a reduced space (205 X 95 cm), both inside and outside of the platform that has been specifically designed to ensure a solid stability on the ground and sized for a convenient use of any bench or step. It also allows a total accessibility even in cases of users with disabilities.



LOAD PROFILES

WITH SINTESI YOU CAN SIMULATE MANY WORKING CONDITIONS, PERSONALIZING YOUR TRAINING, AS NEVER BEFORE



CONSTANT LOAD



ECCENTRIC LOAD



AUXOTONIC LOAD



MAGNETIC LOAD



VISCOUS LOAD



ISOKINETIC LOAD



VIBRANT LOAD



GUIDED MOVEMENT LOAD



ISOMETRIC LOAD

affected by the inertial influence of the load.

AUX - it progressively increases the load depending on the extension of the cable, reproducing the same feeling of an elastic with different levels of resistance.

CST - The load is always constant and is not ECC - It allows you to separate the load in the positive phase from the negative one. You can then increase or decrease the load in the eccentric phase.

VSC - The load varies according to the speed of execution, by setting the viscosity of a fluid. It simulates the same functionality of a rowing machine

ISM - It produces a muscle contraction, even at its top, without any shortening of the muscles and moving of the joint.

VIB - The load varies randomly with adjustable frequency to create instability/ imbalance during the operation.

MAG - The load simulates the same magnetic attraction between two magnets that, the more are close to the zero point, the more the magnetic attraction increases

ISC - The load shows a constant angular velocity motion, regardless of the applied

GMT - It sets the load directly with a divnamic curve according to the extension of the cable, by varying the weight in the proximity of specific points.

STRENGTH HYPERTROPHY



FUNCTIONAL TRAINING

CARDIO



SPORT PERFORMANCE RETURN TO PLAY

PROPRIOCETION STABILITY





TONING SLIMMING

POWER





FUNCTIONAL RECOVER

